

April 2016

## New version of software for Control Omni

Abilia would like to inform you that there will be a new version of the software in Control Omni. It will be released 2016-04-26. The new official version is 1.4.0

New version can be downloaded from Abilia's website 2016-04-26.

<http://www.abilia.com/en/service-and-support/software-updates/software-updates-environmental-control/upgrade-control-omni>

### Version 1.4.0


- Bug fix: Font Size

### Version 1.3.9

- **Enhancement: Full screen scan**  
For each scan step shows each button in full-screen size. The setting is individual for each page and can be set to On or Off.

**Advantage:** Good for visually impaired and you will not be distracted

**Add or remove a Full Screen Scan by using the following procedure:**


 → Pages/buttons → Select page → browse among the pages available →  
Select this page → Properties → Full Screen Scan

- **Enhancement: Automatic start scan**

The scanning will start automatically when coming to a new page. The setting can be set to On or Off and applies to all pages in the user environment.

**Advantage:** Speed up the scanning

**Change Automatic Scanning Start by using the following procedure:**




- **Enhancement: Automatic return**

Returns automatically to Start page if number of scanning turns is completed. The setting can be set to On or Off and applies to all pages in the user environment.

**Advantage:** The scanning will automatically return to the start page if you do not make a choice

**Change Automatic return by using the following procedure:**




- **Enhancement: Font size**

Added ability for large text on button.

**Advantage:** Good for visually impaired and provides cognitive support

**Write text on a button by using the following procedure:**


 → **Pages/buttons** → **Select page** → Select the page on which the button to be edited is located → **Select this page** → **Edit button** → Select the button that is to be edited → **Appearance** → **Change text** or **Change font size**

- **Enhancement: Repetitions of a Macro**

Added more functionality to the Macro function. Now able to repeat a built Macro up to 20 repetitions and the possibility of selecting the delay between repetitions (0-5 sec).

**Advantage:** E.g. if you want to program the channels on the TV to be stepping up automatically.

**Repeat a Macro by using the following procedure:**


 → **Pages/buttons** → **Select page** → Select a page on which the macro is to be stored → **Select this page** → **Edit button** → Select the button on which the macro is to be stored → **Function** → **Environment** → **Macro** → **Start macro recording** → Select the buttons that are to be included in the macro → **Repeat macro**

- **Enhancement: First step delay**

Added "First step delay" as a scanning property. The setting can be set to On or Off.

**Advantage:** No waiting on the first scan step, especially if you have slow scan

**Change First step delay by using the following procedure:**



- **Enhancement: Touch Screen Scan**

The whole touchscreen to become scan input. This can be set to On or Off and applies to all user pages.

**Advantage:** This setting makes the whole touchscreen to become scan input, good when using it in combination of Full Screen Scan

**Change to Touch screen Scan by using the following procedure:**



- **Enhancement: Radio and TV logos**

Introduced UK TV- & Radio-channels, Danish TV channels, and Norwegian TV channels logos pictures.

We recommend to update all Control Omni to this version, not urgent but when possible.

With kind regards

Peter Linghammar

Product manager Environmental Controls

Sergey Saltanov

Export manager

